# LET'S REDUCE SSI TOGETHER

# Applied Medicalの SSIへの取り組み

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“Our meta-analysis found that dual-ring wound protectors reduce the odds of SSI in patients undergoing lower gastrointestinal surgery.”

“We demonstrated evidence of a subgroup difference where dual-ring wound protectors reduced SSIs while single-ring retractors did not, which provides greater insight in the choice of wound protection devices.”

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“Our study suggests that the use of wound protectors decreases the risk of SSI by 45%. Our number needed to treat suggests that only 10 patients would have to be treated intraoperatively with a wound protector to prevent 1 SSI.”

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“In this study the use of barrier wound protection in elective open colorectal resectional surgery resulted in a clinically significant reduction in incisional surgical site infections.”

“There was a significant reduction in the incidence of incisional surgical site infections when the wound protector was used: 3 of 64 (4.7%) vs 15 of 66 (22.7%); P = .004.”

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“Superficial incisional SSI was significantly diminished in the ALEXIS wound retractor group (P=0.006).”

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Wound infection decreased significantly in the With Alexis retractor group by 100% (0% in the With Alexis retractor group vs. 8.1% in the Without Alexis retractor group)

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“[U]se of a plastic wound retractor may result in reduced enteric bacterial colonization of the surgical incision site during gastrointestinal surgery. Reduced colonization of the surgical incision site by enteric bacteria due to the use of a plastic wound retractor should result in a reduction in SSI following gastrointestinal surgery.”

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“Among adult patients with intrabiliary stents, the use of a dual-ring wound protector during pancreaticoduodenectomy significantly reduces the risk of incisional SSI.”

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*The study was terminated early with 109/300 patients after an interim analysis showed a significant 89% reduction in infection when the Alexis retractor was used (1.6% for wound protection arm vs. 14.6% for traditional retraction arm)*

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*"[T]he use of plastic-sheath wound retractors such as the Alexis® O C-Section Retractor compared to the traditional Collins self-retaining metal retractor in low risk women, having the first cesarean is associated with a significantly reduced risk of surgical site infection."*

...

*"There is significant reduction in the use of electric cautery for subcutaneous bleeding, bowel handling and postoperative pain. Operator satisfaction is improved and postoperative pain is less."*

IS ALEXIS  
BRAND  
PART OF YOUR  
STANDARD  
OF CARE?

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© 2022 Applied Medical Resources Corporation. All rights reserved.